



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
High Frequency Mobile X-ray Machine With Battery Back Up

Product Code :
SAMCOMEDICAL430001A010



Description :

High Frequency Mobile X-ray Machine With Battery Back Up Use : NET mobile high-frequency X-ray photography diagnostic equipment can be used in radiology, orthopedics, wards, emergency rooms, operating rooms, and ICU for radiography on human body, such as head, limbs, chest and spine. It has build in battery back up and can work during power supply. It is suitable for radiography requirement in emergency conditions and outdoors. Features : With High voltage generator LCD display, parameter memory, Techniques of kV analog closed loop control, mAs digital closed loop control and Microprocessor based control to ensure precision and repeatability of the dose With 50 preset exposure programs. Process parameters can be changed and modified. Emergency stop button, Screen display of Error number in case of problem Automatic recovery system of parameters on power cut off. With high voltage protection, tube over-current protection, output overload protection function Pillar can rotate 90 deg clockwise and anticlockwise With battery back up. Can take upto 2000 images in case of power failure



Samco Medical

Website: www.samcomedical.com,