



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Ventilator with 10.4

Product Code :
SAMCOMEDICAL4006A001

Description :

Ventilator with 10.4 Screen Work condition : Gas source : O₂ (For medical use) Rated Working Pressure : 280kPa ~ 600kPa Power Source : Voltage : ~100V - 240V, Frequency : 50Hz / 60Hz Ventilation Mode : IPPV, A/C, SIMC, SPONT, CPAP



Samco Medical

Website: www.samcomedical.com,