



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Cervical Collar Hard

Product Code :
SAMCOMEDICAL100003A004

Description :

Cervical Collar Hard

Technical Specification :

Hard Cervical Collar (Adjustable), with velcro locks, sizes small, medium, large. Made of high density polythene sheet. Two Piece height adjustment to maintain desired flexion or extension. Banian Foam Padding for comfort. Velcro closures for easy application, removal and adjustment. Ventilation through holes to allow air circulation for comfort. Washable. Cervical Collar Hard, Cervical Support Exporters, Cervical Support Manufacturers



Samco Medical

Website: www.samcomedical.com,