



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Bilevel Positive Airway Pressure with Multiple Mode

Product Code :
SAMCOMEDICAL4006A005



Description :

Bilevel Positive Airway Pressure with Multiple Mode Application : NET BPAP is mainly used for treating obstructive sleep apnea-hypopnea syndrome (OSAHS) and it is not for patient life support. Features: Work mode : CPAP, S, S/T, T Air leak compensation: Dynamic monitoring of the pressure compensation/ Dynamic monitoring of flow compensation. Large color LCD display: It is very convenient to see the setting parameters and statistical data. It can display waveform of breathing air and pressure in real time. Equipped with Multi-parameter Alarm, Mask off alarm, suffocation alarm, high respiratory rate alarm, low respiratory rate alarm



Samco Medical

Website: www.samcomedical.com,