



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Electrical Anesthetic Ceiling Pendant

Product Code :
SAMCOMEDICAL4007A002



Description :

Electrical Anesthetic Ceiling Pendant Key features : Cross Arm and Control box are molded from High strength Aluminum and the Surface is oxidized for durability. Cross Arm's rotation and height adjustment are controlled by imported electric motor. Rotation is flexible and free. Pneumatic brakes and mechanical brakes provide firm and accurate locking position. Power supply, Gas pipeline and communication lines are arranged separately in one set Equipped with GERMAN brand (GENTEC) standard gas terminal. High weight bearing capacity, Bearing suspension for Anesthesia machine. Anesthesia machine not included. Omni bearing movement makes operation convenient and faster Basic configuration : One Gas discharge terminal for compressed air, oxygen, nitrous oxide, suction and anesthetic gas Four power outlets One network outlets



Samco Medical

Website: www.samcomedical.com,