



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Surgical Ceiling Pendant with Saline Rod Fixation

**Product Code :**  
SAMCOMEDICAL4007A004



**Description :**

Surgical Ceiling Pendant with Saline Rod Fixation Key features : Cross Arm and Control box are molded from high strength Aluminum and the surface is oxidized for durability. Horizontal rotation can accurately position the system. Pneumatic brakes and mechanical brakes provide firm and accurate locking position. Power supply, Gas pipeline and communication line are arranged separately in one set Equipped with GERMAN brand (GENTEC) standard gas terminal Basic configuration : One Gas terminal for compressed air, oxygen and suction Four power outlets One network outlets



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),