



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Surgical Ceiling Pendant with Saline Rod Fixation, Board and Drawer

Product Code :

SAMCOMEDICAL4007A005

**Description :**

Surgical Ceiling Pendant with Saline Rod Fixation, Board and Drawer Key features : Cross Arm and Control box are molded from High strength Aluminum and the Surface is oxidized for durability. Horizontal rotation can accurately position the system. Pneumatic brakes and mechanical brakes provide firm and accurate position locking. Power supply, Gas pipeline and communication line are arranged separately in one set Equipped with GERMAN brand (GENTEC) standard gas terminal Basic configuration : One Gas terminal for compressed air, oxygen and suction Four power outlets One network outlets



Samco Medical

Website: www.samcomedical.com,