



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Surgical Ceiling Pendant with Board and Drawer

Product Code :
SAMCOMEDICAL4007A006



Description :

Surgical Ceiling Pendant with Board and Drawer Key features : Cross Arm and Control box are molded from High strength Aluminum and the Surface is oxidized for durability. Horizontal rotation can accurately position the system. Pneumatic brakes and mechanical brakes provide firm and accurate position locking. Power supply, Gas pipeline and communication line are arranged separately in one set Equipped with GERMAN brand (GENTEC) standard gas terminal. Provides independent area for Monitoring equipments and Infusion liquid. Basic configuration : One Gas terminal for compressed air, oxygen and suction Four power outlets One network outlets



Samco Medical

Website: www.samcomedical.com,