



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Surgical Anesthetic Ceiling Pendant

Product Code :
SAMCOMEDICAL4007A009



Description :

Surgical Anesthetic Ceiling Pendant Key features : Cross Arm and Control box are molded from High strength Aluminum and the Surface is oxidized for durability. Horizontal rotation can accurately position the system. Pneumatic brakes and mechanical brakes provide firm and accurate position locking. Power supply, Gas pipeline and communication line are arranged separately in one set Equipped with GERMAN brand (GENTEC) standard gas terminal. High bearing capacity crane for suspension of Anesthetic trolley Basic configuration : One Gas discharge terminal for compressed air, oxygen, nitrous oxide, suction and anesthetic gas Four power outlets One network outlets



Samco Medical

Website: www.samcomedical.com,