



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Bridge Ceiling Pendant

Product Code :
SAMCOMEDICAL4007A010



Description :

Bridge Ceiling Pendant Key features : Bridge pendant is specially designed for ICU and recovery rooms. It is Equipped with Light in the middle and all kinds of terminal such as medical gas and power network. Free movement makes operation and nursing more convenient. Additional equipment such as communication machine, monitor, view telephone, back ground music, teaching and long distance consultation system can be equipped as per requirement. Cross beam and Control box are molded from High strength Aluminum and the Surface is oxidized for durability. Horizontal rotation is achieved by imported electric motor and can accurately position the system. Pneumatic brakes and mechanical brakes provide firm and accurate position locking. Power supply, Gas pipeline and communication line are arranged separately in one set Equipped with GERMAN brand (GENTEC) standard gas terminal. Provides indePENDANT area for Monitoring equipments and Infusion liquids. Horizontal translation of two hangings provides omni bearing One Gas discharge terminal for compressed air, oxygen, suction. Ten power outlets Two network outlets



Samco Medical

Website: www.samcomedical.com,