



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Medical Column with Two Shelf and Saline Rod

Product Code :
SAMCOMEDICAL4007A011



Description :

Medical Column with Two Shelf and Saline Rod Key features : Economical pendant designed for sick rooms and small operation theatre. Hanging Pole and Control box are molded from High strength Aluminum and the Surface is oxidized for durability. With mechanical brake Power supply, Gas pipeline and communication line are arranged separately in one set Provides two shelf for Monitoring equipments and saline rods for Infusion liquids. Equipped with GERMAN brand (GENTEC) standard gas terminal. Basic configuration : One Gas discharge terminal for compressed air, oxygen, suction Four power outlets One network outlets



Samco Medical

Website: www.samcomedical.com,