



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Medical Column Multifunctional

Product Code :
SAMCOMEDICAL4007A012



Description :

Medical Column Multifunctional Key features : Multifunctional Rotatable pendant designed for sick rooms and small operation theater. This is an ideal workstation for Hanging Pole and Control box are molded from High strength Aluminum and the Surface is oxidized for durability. With mechanical brake Power supply, Gas pipeline and communication line are arranged separately in one set Provides two shelf for Monitoring equipments and saline rods for Infusion liquids. Equipped with GERMAN brand (GENTEC) standard gas terminal. Basic configuration : Two Gas discharge terminal for compressed air, oxygen, suction Four power outlets One network outlets



Samco Medical

Website: www.samcomedical.com,