



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Electrosurgical Unit : 350W

Product Code :
SAMCOMEDICAL4007A013



Description :

Electrosurgical Unit : 350W Standard Matching Accessories: Disposable Electrosurgical Pencil : 5 pcs
Disposable ESU Pad (split) : 5 pcs Disposable ESU Pad (solid) : 5 pcs Cable for Reusable Bipolar Forceps : 1
pc Dual Pedal Switch : 1 pc Individual Pedal Switch : 1 pc Instructions Manual : 1 pc Features : Automatic
Detection Function for Startup and Error Function Memory Function, 10 memory settings, previous power setting
digits available after restart. Neutral electrode contact quality inspection system. Alarm is raised once it detect
that the contact area is getting small to the certain dangerous situation. Large vessel sealing system, With argon
system. With excellent compatibility, it connects with Cystoscopy, arthroscopy, laparoscopy and hysteroscopy
and providing them the power source to cutting and coagulation. It also operated at the urethral resection of the
prostate (TURP), the Arthroscopic surgery and the usage of gasification electrode surgery (TUVF) etc. Each
output terminal is independent output design Transportation Data: Size : 43cm (L) x 38cm (W) x 16.5cm (H)
(UNIT) ; 56 (L) x 63 (W) x 28cm (H) (outer carton) Weight : 7kg (Net Weight) ; 10.5kg (Gross Weight)



Samco Medical

Website: www.samcomedical.com,