



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical Ceiling Shadowless OT Light

Product Code :
SAMCOMEDICAL4007A016



Description :

Samcomedical Ceiling Shadowless OT Light with Single Dome Dia 504mm, Operated by Remote Control

Technical Specification :

Ceiling Shadowless OT Light with Single Dome Dia 504mm, Operated by Remote Control Key Structure A New superior quality O.T. Light executed by a Remote Control. This Light provides the superior performance and is good for all kinds of Surgery. Lights arms are made of Stainless Steel Fitting. The light is equipped with the reserve light source, which is activated automatically in case the main light source brakes down. LED Display that helps the user to adjust the Intensity of Light from 0 to 100% & mode system Technical specifications Dome head : 504mm Halogen Bulbs : 24V x 150W x 2 nos. Light Intensity at 1 mt : 90,000 Lux Height Adjustment : 1078mm Dichoric Glass Reflector : 390mm Color Temperature : 4200 K 300 Temp rise in field : 2 - 3 C Control Unit : SMPS



Samco Medical

Website: www.samcomedical.com,