



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical Air Compressor Nebulizer

Product Code :
SAMCOMEDICAL120001A011



Description :

Samcomedical Air Compressor Nebulizer

Technical Specification :

Compressor Nebulizer delivers aerosol medication therapy effectively and efficiently. This self-contained unit is stylish, attractive. The easy carry handle, two storage compartments and built-in nebulizer holder makes this unit convenient. Power supply: AC220 50HZ Max negative pressure compressing pump: 0.15MPa Free air flux compressing pump: 10L/min Max Nebulizing rate: 0.10mL/min Noise: 60dBA Air Compressor Nebulizer, Emergency Products Exporters, Emergency Products Manufacturers



Samco Medical

Website: www.samcomedical.com,