



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Samcomedical Ceiling Shadowless OT Light Single

Product Code :
SAMCOMEDICAL4007A017



Description :

Samcomedical Ceiling Shadowless OT Light with Single Dome Dia 515mm with 3 Reflectors

Technical Specification :

Ceiling Shadowless OT Light with Single Dome Dia 515mm with 3 Reflectors Key Structure This Light provides the superior performance and is good for all kinds of Surgery. Lights arms are made of Stainless Steel Fitting. Sterilizable handle has been provided in the middle of the dome for regulating the diameter of focus. Light intensity can be adjusted by the Digital Control Panel. Technical specifications Dome head : 515mm Halogen Bulbs : 24V x 70W x 3 nos. Light Intensity at 1 mt : 80,000 Lux Height Adjustment : 1014mm Dichoric Glass Reflector : 200mm Color Temperature : 4200 K 300 Temp rise in field : 2 - 3 C Control Unit : SMPS



Samco Medical

Website: www.samcomedical.com,