



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Ceiling Shadowless OT Light with Single Dome Dia 712mm  
with 7 Reflectors

**Product Code :**

SAMCOMEDICAL4007A019

**Description :**

Ceiling Shadowless OT Light with Single Dome Dia 712mm with 7 Reflectors Key Structure This Operation Reflector Lamp offers the illumination comparable to other standard models. Halogen Lights suitable for Major as well as Minor Surgery, fully spring balanced cool lights. Sterilizable handle has been provided in the middle of the dome to adjust the focus. Parabolic reflectors & color correcting filters. Technical specifications Dome head : 712mm Halogen Bulbs : 24V x 70W x 7 nos. Light Intensity at 1 mt : 90,000 Lux Height Adjustment : 714 mm Parabolic Reflector : 180 mm Color Temperature : 4200 K  $\pm 1/2$  300 Temp rise in field :  $3\frac{1}{2}$  -  $4\frac{1}{2}$  C Control Unit : Transformer



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),