



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Ceiling Shadowless OT Light with Single Dome Dia 550mm  
with 4 Reflectors

**Product Code :**

SAMCOMEDICAL4007A020

**Description :**

Ceiling Shadowless OT Light with Single Dome Dia 550mm with 4 Reflectors Key Structure The light has daylight characters and can be adjusted as per the operation requirements. Lights arms are made of Stainless Steel Fitting. The range of the color temp - 4,200 K  $\pm$  300 as natural color tone. Detachable handle for simplistic sterilization. Parabolic reflectors & color correcting filters. Technical specifications Dome head : 550 mm Halogen Bulbs : 24V x 70W x 4 nos. Light Intensity at 1 mt : 50,000 Lux Height Adjustment : 914 mm Parabolic Reflector : 180 mm Color Temperature : 4200 K  $\pm$  300 Temp rise in field : 2 $\frac{1}{2}$  - 3 $\frac{1}{2}$  C Control Unit : Transformer



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),