



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Ceiling Shadowless OT Light with Single Dome Dia 521mm

Product Code :
SAMCOMEDICAL4007A021



Description :

Ceiling Shadowless OT Light with Single Dome Dia 521mm Key Structure The Light Head provides the superior performance and is good for all kinds of surgeries. Light arms are made of Stainless Steel Fitting. The light is equipped with the reserve light source, which is activated in case the main light source brakes down. Sterilizable handle has been provided on the dome for easy movement. Heat absorbing & color correcting filter glass.
Technical specifications Dome head : 521mm Halogen Bulbs : 24V x 150W x 2 nos. Light Intensity at 1 mt : 50,000 Lux Height Adjustment : 1014mm Color Temperature : 4200 K $\pm \frac{1}{2}$ 300 Temp rise in field : $2\frac{1}{2} - 3\frac{1}{2}$ C
Control Unit : Transformer



Samco Medical

Website: www.samcomedical.com,