



Email : sales@samcomedical.com  
Phone: 1-300-400-8211

**Product Name :**  
Ceiling Shadowless OT Light Single Dome Dia 405mm

**Product Code :**  
SAMCOMEDICAL4007A022



**Description :**

Ceiling Shadowless OT Light Single Dome Dia 405mm Key Structure The light is equipped with reserve light source, which is activated in case the main light source brakes down. Lights arms are made of Stainless Steel Fitting. Sterilizable handle has been provided in the middle of the dome. Heat absorbing & color correcting filter glass. Technical specifications Dome head : 405 mm Halogen Bulbs : 24V x 150W x 2 nos. Light Intensity at 1 mt : 60,000 Lux Height Adjustment : 1014mm Parabolic Reflector : 355mm Color Temperature : 4200 K  $\pm$  300 Temp rise in field :  $2\frac{1}{2}$  -  $3\frac{1}{2}$  C Control Unit : Transformer



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),