



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Pen-Type Digital Thermometer

Product Code :
SAMCOMEDICAL6001A022

Description :

Pen-Type Digital Thermometer

Technical Specification :

Display range: 32C-42C(90F - 107.6F), Accuracy: 0.1 C(0.2F), Min Scale: 0.1, Memory: Last measuring reading, Measurement time: 60½ seconds in oral, 100 seconds underarm, Beep function, Auto shut-off, Battery: 1.5V button battery (LR/SR-41). Pen-Type Digital Thermometer, Diagnostic Equipments Exporters, Diagnostic Equipments Manufacturers



Samco Medical

Website: www.samcomedical.com,