



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Ceiling Shadowless Surgical Operation Light 6 Bulbs - Superior (6 Bulbs, 710mm Dome Dia.)

Product Code :

SAMCOMEDICAL4007A027

**Description :**

Ceiling Shadowless Surgical Operation Light 6 Bulbs - Superior (6 Bulbs, 710mm Dome Dia.) Light head : 710 mm dome dia., 6 bulbs. Light Intensity : 110,000 Lux. Halogen Bulb : 24 V x 50 W x 6 Nos. Power Consumption : 300VA. Action Radius : 1,400 mm. Colour Temperature : 4200 $\pm$ 1/2 K. Temperature Rise : 5 $\pm$ 1/2-8 $\pm$ 1/2C. Average Life of Halogen Bulbs : 4,000 hrs. Power Requirement : 220V AC. Room Height : 9 ft



Samco Medical

Website: www.samcomedical.com,