



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Ceiling Shadowless OT Light with Dia 515mm x 2 Domes,
Each Dome Having 3 Reflectors

Product Code :

SAMCOMEDICAL4007A030

**Description :**

Ceiling Shadowless OT Light with Dia 515mm x 2 Domes, Each Dome Having 3 Reflectors Key Structure This Light provides superior performance and is good for all kinds of Surgery. Inbuilt is an Imported Dichoric Glass Reflector, which increases the Lux up to 80,000 (i.e. 10%). Lights arms are made of stainless Steel Fitting. During the long time operation, there is no fatigue to the eyes, due to the range of the color temp - 4200K i.e. 1/2 300 as natural color tone. Sterilizable handle has been provided in the middle of the dome for regulating the focus. Light intensity can be adjusted by the Digital Control panel. Diffuser spreads the light uniformly. It protects the light and gives excellent shadow reduction. It also reduces infrared radiation and enhances the light's appearance. Technical specifications Dome head : 515 mm x 2 Halogen Bulbs : 24V x 70W x 6 nos. Light Intensity at 1 mt : 1,60,000 Lux Dichoric Glass reflector : 200mm Height Adjustment : 1014mm Color Temperature : 4200 K i.e. 1/2 300 Temp rise in field : 2i.e. 1/2- 3i.e. 1/2 C Control Unit : SMPS



Samco Medical

Website: www.samcomedical.com,