



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Ceiling Shadowless OT Light with Dia 550mm x 2 Domes,  
Each Dome Having 4 Reflectors

**Product Code :**  
SAMCOMEDICAL4007A033



**Description :**

Ceiling Shadowless OT Light with Dia 550mm x 2 Domes, Each Dome Having 4 Reflectors Key Structure The light has daylight characters and can be adjusted as per the operation requirements. Lights arms are made of Stainless Steel Fitting. Detachable handle for simplistic sterilization. Parabolic reflectors & color correcting filter glass. Technical specifications Dome head : 550 mm x 2 Halogen Bulbs : 24V x 70W x 8 nos. Light Intensity at 1 mt : 1,00,000 Lux Height Adjustment : 914mm Parabolic Reflector : 180mm Color Temperature : 4200 K  $\pm 1/2$  300 Temp rise in field :  $2\frac{1}{2}$  -  $3\frac{1}{2}$  C Control Unit : Transformer



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),