



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Ceiling Shadowless OT Light with Dia 521mm x 2 Domes

**Product Code :**  
SAMCOMEDICAL4007A034



**Description :**

Ceiling Shadowless OT Light with Dia 521mm x 2 Domes Key Structure The Light Head provides the superior performance and is good for all kinds of surgeries. The compact light head with shadow reduction control makes this light suitable for wide variety of procedures, including the General Surgery, Micro Surgery, E.N.T., OB/GYN, Orthopedic and Cosmetic Surgery. A Shadowless light suitable for all sizes of Operation Theatres as it gives maximum light output with minimum heat dissipation. Lights arms are made of Stainless Steel Fitting. The range of the color temp - 4,200K  $\pm$  1/2 300 as natural color tone. The light is equipped with the reserve light source, which is activated in case the main light source brakes down. Sterilizable handle has been provided on the dome for easy movement. Technical specifications Dome head : 521mm x 2 Halogen Bulbs : 24V x 150W x 4 nos. Light Intensity at 1 mt : 1,00,000 Lux Height Adjustment : 1014mm Color Temperature : 4200 K  $\pm$  1/2 300 Temp rise in field : 2 $\pm$  1/2 - 3 $\pm$  1/2 C Control Unit : Transformer



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),