



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Ceiling Shadowless Surgical Operation Light Twin Dome -
Superior (4+4 Bulbs)

Product Code :
SAMCOMEDICAL4007A036



Description :

Ceiling Shadowless Surgical Operation Light Twin Dome - Superior (4+4 Bulbs) Light head : 500 mm dome dia., 4 bulbs each. Light Intensity : 75,000 Lux. each. Halogen Bulb : 24V x 50W x 8 Nos. Power Consumption : 400VA. Action Radius : - Upper Arm 1,400 mm. - Lower Arm 1,200 mm. Colour Temperature : 4200 $\pm$ 1/2 K. Temperature Rise : 5 $\pm$ 1/2-8 $\pm$ 1/2C. Average life of halogen bulbs : 4,000 Hrs. Power Requirement : 220V AC. Room Height : 9 ft



Samco Medical

Website: www.samcomedical.com,