



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Ceiling Shadowless OT Light with Dia 504mm X 3 Domes,  
Operated by Remote Control

**Product Code :**

SAMCOMEDICAL4007A038

**Description :**

Ceiling Shadowless OT Light with Dia 504mm X 3 Domes, Operated by Remote Control A new superior quality O.T. Light executed by a wireless Remote Control. This Light provides the superior performance and is good for all kinds of Surgery. Lights arms are made of Stainless Steel Fitting. The light is equipped with the reserve light source, which is activated automatically in case the main light source brakes down. Sterilizable handle is provided in the middle of the dome. Inbuilt is an Imported Dichoric Glass Reflector, which increases the Lux up to 2,70,000 (i.e. 10%). LED Display on the dome that helps the user to adjust the Intensity of Light from 0 to 100% & mode system. Technical specifications Dome head : 504 mm x 3 Halogen Bulbs : 24V x 150W x 6 nos. Light Intensity at 1 mt : 2,70,000 Lux Height Adjustment : 978 mm Dichoric Glass reflector : 390 mm Color Temperature : 4200 K i.e. 300 Temp rise in field (1 mt) : 2 i.e. 3 i.e. C Control Unit : SMPS



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),