



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Ceiling Shadowless OT Light with Dia 405mm x 3 Domes

Product Code :
SAMCOMEDICAL4007A040



Description :

Ceiling Shadowless OT Light with Dia 405mm x 3 Domes This Light provides the superior performance and is good for all kinds of Surgery. The light is equipped with the reserve light source, which is activated in case the main light source brakes down. Lights arms are made of Stainless steel Fitting. Sterilizable handle has been provided in the middle of the dome. Inbuilt parabolic reflector, which increases the Lux up to 1,80,000 (i.e. 10%). Heat absorbing & color correcting filter glass. Technical specifications Dome head : 405 mm x 3 Halogen Bulbs : 24V x 150W x 6 nos. Light Intensity at 1 mt : 1,80,000 Lux Height Adjustment : 1014mm Parabolic Reflector : 355 mm Color Temperature : 4200 K i.e. 1/2 300 Temp rise in field : 2 i.e. 1/2- 3 i.e. 1/2 C Control Unit : Transformer



Samco Medical

Website: www.samcomedical.com,