



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Wall Mounted Shadowless OT Lights with Single Dome Dia 465mm

**Product Code :**

SAMCOMEDICAL4007A041

**Description :**

Wall Mounted Shadowless OT Light with Single Dome Dia 465mm Key Structure Light head is equipped with two bulbs (the primary bulb and the back-up bulb). When the primary bulb fails, the back-up bulb lights up automatically to ensure continuous uninterrupted lighting throughout the case. Lights arms are made of Stainless Steel Fitting. Sterilizable handle has been provided in the middle of the dome in feather touch moment. Digital Control Panel to increase or decrease intensity of the light. Technical specifications Dome head : 465mm Halogen Bulbs : 24V x 150W x 2 nos. Light Intensity at 1 mt : 80,000 Lux Height Adjustment : 1065mm Dichoric Glass Reflector : 365mm Color Temperature : 4200 K  $\pm 1/2$  300 Temp rise in field :  $2\frac{1}{2} - 3\frac{1}{2}$  C Control Unit : SMPS



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),