



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Wall Mounted Shadowless OT Light with Single Dome Dia 515mm with 3 Reflectors

Product Code :

SAMCOMEDICAL4007A042

**Description :**

Wall Mounted Shadowless OT Light with Single Dome Dia 515mm with 3 Reflectors Key Structure This Light provides the superior performance and is good for all kinds of Surgery. Lights arms are made of Stainless Steel Fitting. Sterilizable handle has been provided in the middle of the dome for regulating the diameter of focus. Light intensity can be adjusted by the Digital Control Panel. Diffuser spreads the light uniformly. It protects the light and gives excellent shadow reduction. It also reduces infrared radiation and enhances the light's appearance. Technical specifications Dome head : 515mm Halogen Bulbs : 24V x 70W x 3 nos. Light Intensity at 1 mt : 80,000 Lux Height Adjustment : 1014mm Dichoric Glass Reflector : 200mm Color Temperature : 4200 K $\frac{1}{2}$ 300 Temp rise in field : $2\frac{1}{2}$ - $3\frac{1}{2}$ C Control Unit : SMPS



Samco Medical

Website: www.samcomedical.com,