



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**

Mobile Shadowless OT Light with Single Dome Dia 504mm,  
Operated by Remote Control

**Product Code :**

SAMCOMEDICAL4007A044

**Description :**

Mobile Shadowless OT Light with Single Dome Dia 504mm, Operated by Remote Control Key Structure A New superior quality O.T. Light executed by a Remote Control. This Light provides the superior performance and is good for all kinds of Surgery. Light arms are made of Stainless Steel Fitting. The light is equipped with the reserve light source, which is activated automatically in case the main light source brakes down. Sterilizable handle has been provided in the middle of the dome. LED Display that helps the user to adjust the Intensity of Light from 0 to 100% & mode system. Technical specifications Dome head : 504mm Halogen Bulbs : 24V x 150W x 2 nos. Light Intensity at 1 mt : 90,000 Lux Height Adjustment : 1078mm Dichoric Glass Reflector : 390mm Color Temperature : 4200 K  $\pm \frac{1}{2}$  300 Temp. rise in field :  $2\frac{1}{2}$  -  $3\frac{1}{2}$  C Control Unit : SMPS



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),