



Email : sales@samcomedical.com  
Phone: 1-300-400-8211

**Product Name :**

Mobile Shadowless OT Light with Single Dome Dia 465mm

**Product Code :**

SAMCOMEDICAL4007A046

**Description :**

Mobile Shadowless OT Light with Single Dome Dia 465mm Light head is equipped with two bulbs (the primary bulb and the back-up bulb). When the primary bulb fails, the back-up bulb lights up automatically to ensure continuous uninterrupted lighting throughout the case. Light arms are made of Stainless Steel Fitting. Sterilizable handle has been provided in the middle of the dome in feather touch moment. Digital Control Panel to increase or decrease the intensity of the light. Technical specifications Dome head : 465mm Halogen Bulbs : 24V x 150W x 2 nos. Light Intensity at 1 mt : 80,000 Lux Height Adjustment : 1065mm Dichoric Glass Reflector : 365mm Color Temperature : 4200 K  $\pm 1/2$  300 Temp. rise in field :  $2\frac{1}{2}$  -  $3\frac{1}{2}$  C Control Unit : SMPS



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),