



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Baby Nipple Thermometer

Product Code :
SAMCOMEDICAL130001A014

Description :

Baby Nipple Thermometer

Technical Specification :

Baby Nipple Thermometer - Display range: 34C~44.0C(93.2F~111.2F, Accuracy: .1C(0.2F), Min Scale: 0.1, Fast,accurate and easy operation, Beep function, Digital display, Baby safe, Mercury free, Battery: 1.5V button battery(LR/SR41), Battery life: about 100 hours under continuous operation. Baby Nipple Thermometer, Baby Care Equipments Exporters, Baby Care Equipments Manufacturers



Samco Medical

Website: www.samcomedical.com,