



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Mobile Shadowless OT Light with Single Dome dia 550mm
with 4 Reflectors

Product Code :
SAMCOMEDICAL4007A047



Description :

Mobile Shadowless OT Light with Single Dome dia 550mm with 4 Reflectors Key Structure The light has daylight characters and can be adjusted as per the operation requirements. Light arms are made of Stainless Steel Fitting. The range of the color temp. - 4200 K \pm 300 as natural color tone. Detachable handle for simplistic sterilization. Parabolic reflectors & color correcting filters. Technical specifications Dome head : 550mm Halogen Bulbs : 24V x 70W x 4 nos. Light Intensity at 1 mt : 50,000 Lux Height Adjustment : 914mm Parabolic Reflector : 180mm Color Temperature : 4200 K \pm 300 Temp. rise in field : 2 $^{\circ}$ \pm 1/2 - 3 $^{\circ}$ \pm 1/2 C Control Unit : Transformer



Samco Medical

Website: www.samcomedical.com,