



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Mobile Shadowless OT Light with Single Dome Dia 521mm

Product Code :

SAMCOMEDICAL4007A048



Description :

Mobile Shadowless OT Light with Single Dome Dia 521mm Key Structure The Light Head provides superior performance and is good for all kinds of surgeries. The compact light head with shadow reduction control makes this light suitable for wide variety of procedures, including General Surgery, Micro Surgery, E.N.T., OB/GYN, Orthopedic and Cosmetic Surgery. A Shadow less light suitable for all sizes of Operation Theatres as it gives maximum light output with minimum heat dissipation. Light arms are made of Stainless Steel Fitting. The range of the color temp. - 4,200 K \pm 300 as natural color tone. The light equipped with reserve light source, which is activated in case the main light source brakes down. Sterilizable handle has been provided on the dome for easy movement. Technical specifications Dome head : 521mm Halogen Bulbs : 24V x 150W x 2 nos. Light Intensity at 1 mt : 50,000 Lux Height Adjustment : 1014mm Color Temperature : 4200 K \pm 300 Temp. rise in field : 2 $\frac{1}{2}$ -3 $\frac{1}{2}$ C Control Unit : Transformer



Samco Medical

Website: www.samcomedical.com,