



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Shadowless Mobile Surgical Operation Light Single Dome - Superior (480mm dome dia.)

Product Code :

SAMCOMEDICAL4007A051

**Description :**

Shadowless Mobile Surgical Operation Light Single Dome - Superior (480mm dome dia.) The mobile surgical light is an improved lighting system designed to deliver cool and sufficient light with a single unit. Most suitable for minor surgical operations and private clinics. Totally enclosed, mirror-finish metallised dome of 480 mm dia. Fitted with heat absorbing, colour correcting cylindrical glass filter and 24V halogen bulb to ensure glare free, cool and brilliant light with optimized colour rendering. The spring balance system allows unprecedented and effortless lighthouse movements. Special halogen bulb (24V x 150W) ensures higher illumination (30,000 Lux) and lesser power consumption in comparison to incandescent bulb. A handy focusing arrangement for proper intensity of light. Provided with a new generation R-CORE step down transformer working on mains 220V AC for 24V supply. Mounted on a heavy, tilt proof steel base with four double wheel castors for silent and smooth floor movements. Beautifully finished with extra hard oven baked polychromatic paint



Samco Medical

Website: www.samcomedical.com,