



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Vision Drum

Product Code :
SAMCOMEDICAL240001A007



Description :

Test: Test chart contain six test Dots, English, Any regional language, Numbers or Pictures, 'C or 'E' Chart Spot light For Retinoscopy, Colour Test, Duochrome Test. Operation : All Functions can be operated by cord-less remote control, i.e. Chart Up Chart Down Vision Test light, Spot light, Colour Test light Duochrome Test light. Cord-less remote works from any direction. Remote : Drum is supplied with Two remote controls. One Cord-less and One with Cord. Cord remote Available for stand by and emergency use. When cord-less remote is not working properly. If you require or need you can use both the remote simultaneously. Mounting : Suitable mounting arrangements for Ophthalmic Chair Unit Mounting and Wall Mounting. Vision Drum, Ophthalmic Equipments Exporters, Ophthalmic Equipments Manufacturers



Samco Medical

Website: www.samcomedical.com,