



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Schiotz Tonometer

Product Code :
SAMCOMEDICAL240001A014



Description :

Schiotz Tonometer: Highly precise eye tonometer for measuring intra-ocular pressure. Its reliability and durability makes this quality measuring instrument indispensable. High quality agate bearing for an extremely long service life. Precision measurement on a scale of 0 to 20 subdivisions and 0 to -1 sub-divisions. Perfect reading of the scale with red pointer. All vital parts are made of stainless steel, other parts are chrome-plated. Complete in a case with velvet-look inserts. Supplied with three weights (5.5 g, 7.5 g, 10 g) and a conversion table. Schiotz Tonometer, Ophthalmic Equipments Exporters, Ophthalmic Equipments Manufacturers



Samco Medical

Website: www.samcomedical.com,