



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Slit Lamp Zoom

Product Code :
SAMCOMEDICAL240001A016



Description :

Zoom Slit Optical Head with Beam Splitter 'C' Mount Slit lamp Zoom Mag with CCD Digital Camera. It also has the features of slit length & width adjustment. This is helpful in all kind of intensity control of illumination i.e low, medium & high. View Type : Gem Zoom Optical Head with Beam Splitter 'C' Mount & Camera Magnification: 5x,10x,15x,20x & 25x continuous zoom Interpupillary Distance : 52mm to 80mm Camera : CCD Digital Camera Illumination Section Slit Width Adjustment : 0 to 12mm without steps Slit Length Adjustment: 0 to 12mm continuous Diameter of Diaphragm :12mm, 10mm, 8mm, 5mm, 3mm, 2mm, 1mm and 0.2mm Angle of Rotation of Slit: 90 Degree Tilt of Slit (Degree): Horizontal in steps of 0, 5, 10, 15 Vertical in steps (Degree) 0, 5, 10, 15, 20 Filters : Cobalt Blue, Green, Neutral Density, Heat absorbing Light Source: Halogen Lamp 12v-50w Intensity Control of Illumination: Low, Medium and High Slit Lamp Zoom, Ophthalmic Equipments Exporters, Ophthalmic Equipments Manufacturers



Samco Medical

Website: www.samcomedical.com,