



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Wheelchair Folding Attendant Type

Product Code :
SAMCOMEDICAL10027A001

Description :

Wheelchair Folding Attendant Type

Technical Specification :

Wheel chair folding, fix armrest pushing (Attendant) type. 12.5cm castor in front & 31cm wheel at the rear with brakes. This wheel chair is designed for general purpose use at Hospitals/ Nursing homes/ Railways/Airports/Shopping malls. Wheelchair Folding Attendant Type, Wheel Chairs Exporters, Wheel Chairs Manufacturers



Samco Medical

Website: www.samcomedical.com,