



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Wheelchair (Economy) Model

**Product Code :**  
SAMCOMEDICAL10027A002

**Description :**

Wheelchair (Economy) Model

**Technical Specification :**

Wheelchair (Economy) Model Folding Model. Made by E.R.W. Tubes 22.22 mm 18 G. Seat & Back - Cushion  
Seat & Back. Wheel Diameter : Rear 24 x 1.5", (without brack, Front Wheel - Castors wheel 4in, Foot Rest -  
E.R.W. Tubes / M.S. Sheet. Wheelchair (Economy) Model, Wheel Chairs Exporters, Wheel Chairs  
Manufacturers Colour : Black/white/Gray as per requirement.



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),