



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Wheelchair Folding (Travel)

Product Code :
SAMCOMEDICAL10027A005

Description :

Wheelchair Folding (Travel)

Technical Specification :

Wheel chair folding, fix armrest, swing away footrest (pushing type) can easily be stored & can be carried in small car. 12.5cm castors in front & 31cm wheel at rear with brakes. This economical model is specially designed for traveling purpose. Wheelchair Folding (Travel), Wheel Chairs Exporters, Wheel Chairs Manufacturers



Samco Medical

Website: www.samcomedical.com,