



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Wheelchair Folding (Travel) Standard

Product Code :
SAMCOMEDICAL10027A006

Description :

Wheelchair Folding (Travel) Standard

Technical Specification :

Wheel chair folding, detachable armrest 15cm castor in front, 31cm wheels at rear with brakes, swing away foot rest for close approach & compact to fit in small car. This standard wheel chair specially design for out door traveling purpose. Wheelchair Folding (Travel) Standard, Wheel Chairs Exporters, Wheel Chairs Manufacturers



Samco Medical

Website: www.samcomedical.com,