



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Wheelchair Folding (Travel) Superior

**Product Code :**  
SAMCOMEDICAL10027A007

**Description :**

Wheelchair Folding (Travel) Superior

**Technical Specification :**

Wheel chair folding detachable armrest. 20cm castors in front with ball bearing movement & 41cm wheels with grey tyre at rear with brakes, swing away and adjustable foot rest for close approach can easily be stored and carried in small car. This superior compact & lightweight chair is suitable for traveling purpose & Rough surface. Wheelchair Folding (Travel) Superior, Wheel Chairs Exporters, Wheel Chairs Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),