



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Wheelchair Non-Folding With Headrest

Product Code :

SAMCOMEDICAL10027A009

Description :

Wheelchair Non-Folding With Headrest

Technical Specification :

Wheelchair Non-Folding With Headrest - Wheel chair Non-folding 15cm castors in front, 31cm wheels at rear with brakes. Reclining back with anatomic headrest. Swing away and elevating detachable leg rest. Wheelchair Non-Folding With Headrest, Wheel Chairs Exporters, Wheel Chairs Manufacturers



Samco Medical

Website: www.samcomedical.com,