



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :

Wheelchair Folding with Reclining Back

Product Code :

SAMCOMEDICAL10027A010

Description :

Wheelchair Folding with Reclining Back

Technical Specification :

Wheel chair folding. Long adjustable reclining back (also provisional head rest on demand) 20cm castor in front, 41cm at rear with brakes, detachable long armrest with side sport. Swing away and Elevating adjustable leg rest. This bed type wheel chair is also suitable for small car Wheelchair Folding with Reclining Back, Wheel Chairs Exporters, Wheel Chairs Manufacturers



Samco Medical

Website: www.samcomedical.com,