



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Hospital Bed With Built-In Backrest

Product Code :
SAMCOMEDICAL1001A002

Description :

Hospital Bed With Built-In Backrest

Technical Specification :

Size : 198 L x 90 W cm x 60 H cms - Frame work made of strong rectangular steep tubes. Two section of CRCA M.S. Sheet uniformly perforated. Backrest section adjustable on ratchet mechanism. Head and leg bows of steel tubular pipes assembled to laminated pannels. Legs of ERW tubular peipes mounted on four heavy duty rubber shoes. Pre-treated and epoxy powder coated.

Optional Accessories : * Side Railing * Saline rod * Mosquito Net Poles

Note : Above Bed can also supplied with castors.



Samco Medical

Website: www.samcomedical.com,