



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Wheelchair Folding with Pot (Steel)

Product Code :
SAMCOMEDICAL10027A012

Description :

Wheelchair Folding with Pot (Steel)

Technical Specification :

Wheel chair folding with pot arrangement detachable arm rest 12.5cm castor in front, 15cm wheel at rear with brakes. This chair is made up of Stainless Steel and is specially designed for bathing purpose (Overall size 23½" x 12½" x 18½"). This model goes directly on the W.C. Wheelchair Folding with Pot (Steel), Wheel Chairs Exporters, Wheel Chairs Manufacturers



Samco Medical

Website: www.samcomedical.com,