



Email : [sales@samcomedical.com](mailto:sales@samcomedical.com)  
Phone: 1-300-400-8211

**Product Name :**  
Wheelchair Self Propelling

**Product Code :**  
SAMCOMEDICAL10027A014

**Description :**

Wheelchair Self Propelling

**Technical Specification :**

Wheelchair Self Propelling - Folding with detachable arm rest ,swing away & detachable adjustable footrest. 20cm castor in front ,61cm wheel at rear. This wheel is design to set the position of legs at any angle. This model is also available in mag wheel. Wheelchair Self Propelling, Wheel Chairs Exporters, Wheel Chairs Manufacturers



**Samco Medical**

Website: [www.samcomedical.com](http://www.samcomedical.com),