



Email : sales@samcomedical.com
Phone: 1-300-400-8211

Product Name :
Wheelchair with Lifting Arrangement

Product Code :
SAMCOMEDICAL10027A015

Description :

Wheelchair with Lifting Arrangement

Technical Specification :

Wheelchair with Lifting Arrangement - Wheel chair folding with flip up armrest, narrow design. Overall width 18 inch (approx) with lifting arrangement this wheel chair is specially designed for pass through narrow passages that are in aircrafts railways. Wheelchair with Lifting Arrangement, Wheel Chairs Exporters, Wheel Chairs Manufacturers



Samco Medical

Website: www.samcomedical.com,